

Intensive Care Unit

Inpatient Care

Neonatology

Home Care

Transafix Fastening Strip No. 1

Often times, gastric or nasal probes as well as leads or drainages are disturbing the patient comfort due to the way they are placed. With our Transafix Fastening Strip No. 1 we developed a solution so the probe or lead doesn't hang loosely on the body.


The Transafix Fastening Strip No. 1 is the perfect solution for fixing smaller volumes, also for children. The basic holder for the Velcro strap, with which the probe or leads are fixed securely, is the Transpore carrier patch.



The lead or probe can either move freely within the Velcro strap or be fixed securely by taking off the protective foil, uncovering an adhesive film. The Velcro strap can be opened and closed as often as you like. Hence, the fastening strip doesn't have to be removed when the probe is exchanged.

The adhesive is breathable and gentle to the skin and can stay on the skin for a longer period of time, even on sensitive skin. It is adaptable and has a strong hold but can also be removed without any residue.

Our Transafix Fastening strip is available in boxes of 50 pieces (sterile) or 100 pieces (non sterile) each.



Benefits

- Can move freely or be fixed in the Velcro strap
- Gentle even to sensitive skin
- Suitable for children and adults
- Perfect for smaller volumes
- Suitable for a long-term treatment
- Maximum patient comfort
- Opening and closing of the Velcro strap possible at any time

Transafix Fastening Strip No. 1

Sterile, with Velcro strap

Article-No. 8800100
PPN*: 111126903089

Non sterile, with Velcro strap

Article-No. 8700150
PPN*: 110050611354

* Pharmacy Product Number

Transafix Fastening Strip No. 1

Overview

Intensive Care Unit

Inpatient Care

Neonatology

Home Care

- Optimal for smaller volumes
- Suitable for long-term treatments
- Fastening strip does NOT have to be removed when exchanging the probe due to the reclosable Velcro strap

Field of application:

Probes, leads, drainages

Application duration:

5-7 days

Illustration




Illustration of the product in original size

**Transafix Fastening Strip No. 1****Sterile, with Velcro strap**

Article-No. 8800100

PPN*: 111126903089

Transafix Fastening Strip No. 1**Non sterile, with Velcro strap**

Article-No. 8700150

PPN*: 110050611354

* Pharmacy Product Number