

Fixing of Endotracheal Tubes

Intensive Care Unit

Emergency Service

Especially when using tubes, a secure and comfortable fixation is important. With our Transafix Tube Fixing System, we developed a perfect solution for this issue.

The set consists of two parts, the 2 cheek plates, available in fleece or hydrocolloid, as well as the fixing strip with a length of 40 cm. The strip can be adjusted individually and is therefore suitable for all different shapes and sizes.

The cheek plates and the fixing strip are connected via a secure Velcro strap. Due to the skin sensitive material, the set can be used for up to 12 days.

Both materials, fleece and hydrocolloid are very gentle to the skin. For patients with sensitive or damaged skin we recommend the use of the hydrocolloid cheek plates.

The cheek plates and fixing strips are available as a set (2 cheek plates and 1 fixing strip) or as individual pieces.



Transafix Tube Fixing System - hydrocolloid

Sterile, Set of 2 cheek plates and 1 fixing strip

Article-No.	8600610
PPN*	110948200634

Transafix Tube Fixing System - fleece

Set of 2 cheek plates and 1 fixing strip

sterile	Article-No.	8800600
non sterile	Article-No.	8700600
	PPN*	110050609994

2 Cheek plates - fleece

sterile	Article-No.	8800700
non sterile	Article-No.	8700700
	PPN*	110050613610

1 Fixing strip

sterile	Article-No.	8800800
	PPN*	110948201200
non sterile	Article-No.	8700800
	PPN*	110050608207

* Pharmacy Product Number

Benefits

- Highest quality
- Easy handling
- Gentle, even to sensitive skin
- Quick and easy cleaning of the mouth possible
- Fixing possible without lifting head
- More economical than common solutions
- No residue
- Application duration of up to 12 days

Fixing of Endotracheal Tubes

Overview

Intensive Care Unit

Emergency Service

- Easy and secure fixing due to the Velcro strap on the cheek plates and the fixing strip
- The fixing can be performed without having to lift the head
- Fits children and adults as well as all sizes of tubing

Field of application:

Endotracheal tubes, nasal probes

Application duration:

up to 12 days

Illustration







Illustration of the products in original size



Transfix cheek plate - hydrocolloid



Transfix cheek plate - fleece



Transfix fixing strip