

Intensive Care Unit

Inpatient Care

Home Care

Transafix Fastening Strip No. 3

Bladder catheters are used in all different kinds of diseases. A discreet use as well as an unrestricted mobility of the patient is necessary for medical and practical reasons. Therefore, it is especially significant to be able to fix the catheter so it is secured on the body but also adjustable to the movement of the patient. Our Transafix Fastening Strip No. 3 offers the perfect solution.




The basic carrier of the Velcro strap is the Transpore carrier patch. It contains an integrated Velcro Strap with which the medical tubing can be fixed securely.

The catheter can either move freely within the Velcro strap or be fixed securely by taking off the protective foil, uncovering an adhesive film. The Velcro strap can be opened and closed as often as you like. Hence, the fastening strip doesn't have to be removed when exchanging the catheter.

The adhesive is breathable and gentle to the skin and can stay on the skin for a longer period of time, even on sensitive skin. It is adaptable and has a strong hold but can also be removed without any residue.

Transafix Fastening No. 3 is available in boxes with 30 pieces (sterile) or 25 pieces (non-sterile) each.



Benefits

- Can move freely or be fixed in the Velcro strap
- Gentle even to sensitive skin
- Perfect for catheters
- Suitable for a long-term treatment
- Maximum patient comfort
- Highest quality
- Opening and closing of the Velcro strap possible at any time

Transafix Fastening Strip No. 3

Sterile, with Velcro strap | 15 x 4,5 cm

Article-No. 8800300

PPN*: 111126909954

Non sterile, with Velcro strap | 15 x 4,5 cm

Article-No. 8700350

PPN*: 111126896442

* Pharmacy Product Number

Transafix Fastening Strip No. 3

Intensive Care Unit

Inpatient Care

Home Care

Overview

- Perfect for larger volumes such as catheters
- Suitable for long-term treatments
- Fastening strip does NOT have to be removed when exchanging the probe due to the resclosable Velcro strap

Field of application:

Bladder catheters, suprapubic catheter as well as wound drainages, central venous catheters, leads, etc.

Application duration:

5-7 days

Illustration





Illustration of the product in original size

**Transafix Fastening Strip No. 3**

Sterile, with Velcro strap | 15 x 4,5 cm

Article-No. 8800300

PPN*: 111126909954

Transafix Fastening Strip No. 3

Non sterile, with Velcro strap | 15 x 4,5 cm

Article-No. 8700350

PPN*: 111126896442

* Pharmacy Product Number