

Venus



Vaginal rejuvenation &
pelvic floor treatments

Limmer
L A S E R 



One handpiece, many configurations



- **Extremely fast treatment** (≥ 0.1 ms / spot)
- **Stacking** of up to 10 levels per fractional dot
- **Flexible dot spacing** (e.g. 200 - 2000 μ m, opt. randomization)
- **Adjustable pattern sizes** (e.g. 1 - 12 mm)
- **Fractional CombiScan** for superficial and deep dots in one shot

Easily to exchange 360° and 45° mirror tips



Intravaginal & external skin treatments with one single scanner system

Typical applications

- Vaginal atrophy
- Stress incontinence
- Postpartum perineal trauma
- Vaginal rejuvenation
- Skin treatments



Scanner remote control with four handpiece buttons



Extractable lens holder for easy cleaning and sterilization

Efficient treatments, easy to handle



Air flow inlet keeps the lens in optimal conditions

Medical degree stainless steel, easy to clean & sterilize

Large scanner patterns for rapid treatments